



TMEM2 Antibody, Biotin conjugated

Product Code	CSB-PA23279D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHN6
Immunogen	Recombinant Human Cell surface hyaluronidase protein (734-866AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface hyaluronidase (EC 3.2.1.35) (Transmembrane protein 2), TMEM2, KIAA1412
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.