



TMEM2 Antibody, Biotin conjugated

Product Code	CSB-PA23279D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHN6
Immunogen	Recombinant Human Cell surface hyaluronidase protein (734-866AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface hyaluronidase (EC 3.2.1.35) (Transmembrane protein 2), TMEM2, KIAA1412
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM2