



ATP2B4 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA23389C0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P23634 |
| Immunogen | Recombinant Human Plasma membrane calcium-transporting ATPase 4 protein (281-345AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Plasma membrane calcium-transporting ATPase 4 (PMCA4) (EC 3.6.3.8) (Matrix-remodeling-associated protein 1) (Plasma membrane calcium ATPase isoform 4) (Plasma membrane calcium pump isoform 4), ATP2B4, ATP2B2, MXRA1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | ATP2B4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |