




tp63 Antibody

Product Code	CSB-PA23919A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E9QG65
Immunogen	Recombinant Zebrafish Tumor protein p63 protein (41-153AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tp63
Immunogen Species	Zebrafish
Research Area	Others
Target Names	tp63

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: tp63 antibody at 7.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 43 kDa

Observed band size: 43, 52 kDa