




rbm47 Antibody

Product Code	CSB-PA23979A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2BIC5
Immunogen	Recombinant Zebrafish RNA-binding motif protein 47 protein (390-524AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	rbm47
Immunogen Species	Zebrafish
Research Area	Others
Target Names	rbm47

Image



Western Blot

Positive WB detected in: Zebrafish tissue, Zebrafish eye tissue

All lanes: rbm47 at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 66 kDa

Observed band size: 66 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.