




mturn Antibody

Product Code	CSB-PA24389A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C8M4
Immunogen	Recombinant Zebrafish Maturin protein (1-133AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Maturin (Maturin neural progenitor differentiation regulator protein homolog), mturn
Immunogen Species	Zebrafish
Research Area	Others
Target Names	mturn

Image



Western Blot

Positive WB detected in: Zebrafish eye tissue

All lanes: mturn at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 16 kDa

Observed band size: 16 kDa