



mturn Antibody, FITC conjugated

Product Code	CSB-PA24389C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C8M4
Immunogen	Recombinant Zebrafish Maturin protein (1-133AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Maturin (Maturin neural progenitor differentiation regulator protein homolog), mturn
Immunogen Species	Zebrafish
Research Area	Others
Target Names	mturn
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.