



th Antibody, Biotin conjugated

Product Code	CSB-PA24419D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9W7K7
Immunogen	Recombinant Zebrafish Tyrosine hydroxylase protein (1-110AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	th
Immunogen Species	Zebrafish
Research Area	Others
Target Names	th
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.