



GLDC Antibody, Biotin conjugated

Product Code	CSB-PA328400LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23378
Immunogen	Recombinant Human Glycine dehydrogenase (decarboxylating), mitochondrial protein (868-984AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycine dehydrogenase (decarboxylating), mitochondrial (EC 1.4.4.2) (Glycine cleavage system P protein) (Glycine decarboxylase) (Glycine dehydrogenase (aminomethyl-transferring)), GLDC, GCSP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GLDC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.