



NS Antibody, Biotin conjugated

Product Code	CSB-PA356065LD01IIE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03495
Immunogen	Recombinant Influenza A virus Non-structural protein 1 protein (152-237AA)
Raised In	Rabbit
Species Reactivity	Influenza A virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-structural protein 1 (NS1) (NS1A), NS
Immunogen Species	Influenza A virus (strain A/Udorn/307/1972 H3N2)
Research Area	Others
Target Names	NS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.