



MOCS2 Antibody, HRP conjugated

Product Code	CSB-PA529746LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O96033
Immunogen	Recombinant Human Molybdopterin synthase sulfur carrier subunit protein (1-88AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Molybdopterin synthase sulfur carrier subunit (MOCO1-A) (Molybdenum cofactor synthesis protein 2 small subunit) (Molybdenum cofactor synthesis protein 2A) (MOCS2A) (Molybdopterin-synthase small subunit) (Sulfur carrier protein MOCS2A), MOCS2, MOCO1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MOCS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.