



ACHE Antibody, FITC conjugated

Product Code	CSB-PA5619C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22303
Immunogen	Recombinant Human Acetylcholinesterase protein (53-162AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acetylcholinesterase (AChE) (EC 3.1.1.7), ACHE
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ACHE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.