



MGAT5 Antibody, Biotin conjugated

Product Code	CSB-PA600992LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q09328
Immunogen	Recombinant Human Alpha-1,6-mannosylglycoprotein 6-beta-N-acetylglucosaminyltransferase A protein (107-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-1,6-mannosylglycoprotein 6-beta-N-acetylglucosaminyltransferase A (EC 2.4.1.155) (Alpha-mannoside beta-1,6-N-acetylglucosaminyltransferase) (GlcNAc-T V) (GNT-V) (Mannoside acetylglucosaminyltransferase 5) (N-acetylglucosaminyl-transferase V), MGAT5, GGNT5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MGAT5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.