



# NEUROD2 Antibody, Biotin conjugated

Product Code	CSB-PA614900LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15784
Immunogen	Recombinant Human Neurogenic differentiation factor 2 protein (254-339AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic differentiation factor 2 (NeuroD2) (Class A basic helix-loop-helix protein 1) (bHLHa1) (NeuroD-related factor) (NDRF), NEUROD2, BHLHA1 NDRF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NEUROD2