



NPRL3 Antibody, Biotin conjugated

Product Code	CSB-PA619643LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12980
Immunogen	Recombinant Human GATOR complex protein NPRL3 protein (349-482AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GATOR complex protein NPRL3 (-14 gene protein) (Alpha-globin regulatory element-containing gene protein) (Nitrogen permease regulator 3-like protein) (Protein CGTHBA), NPRL3, C16orf35 CGTHBA MARE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NPRL3