



NPRL3 Antibody, Biotin conjugated

Product Code	CSB-PA619643LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12980
Immunogen	Recombinant Human GATOR complex protein NPRL3 protein (349-482AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GATOR complex protein NPRL3 (-14 gene protein) (Alpha-globin regulatory element-containing gene protein) (Nitrogen permease regulator 3-like protein) (Protein CGTHBA), NPRL3, C16orf35 CGTHBA MARE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NPRL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.