



CAPS Antibody, HRP conjugated

Product Code	CSB-PA622689LB01HU
Abbreviation	Calcyphosin
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13938
Immunogen	Recombinant Human Calcyphosin protein (1-189AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Calcium-binding protein. May play a role in cellular signaling events (Potential).
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcyphosin (Calcyphosine), CAPS
Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAPS