



CCDC6 Antibody, Biotin conjugated

Product Code	CSB-PA623068LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16204
Immunogen	Recombinant Human Coiled-coil domain-containing protein 6 protein (347-440AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil domain-containing protein 6 (Papillary thyroid carcinoma-encoded protein) (Protein H4), CCDC6, D10S170 TST1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CCDC6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.