



HIC1 Antibody, HRP conjugated

Product Code	CSB-PA623921LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14526
Immunogen	Recombinant Human Hypermethylated in cancer 1 protein (297-418AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hypermethylated in cancer 1 protein (Hic-1) (Zinc finger and BTB domain-containing protein 29), HIC1, ZBTB29
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.