



HIC1 Antibody, Biotin conjugated

Product Code	CSB-PA623921LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14526
Immunogen	Recombinant Human Hypermethylated in cancer 1 protein (297-418AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hypermethylated in cancer 1 protein (Hic-1) (Zinc finger and BTB domain-containing protein 29), HIC1, ZBTB29
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIC1