



ICE2 Antibody, Biotin conjugated

Product Code	CSB-PA717659LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q659A1
Immunogen	Recombinant Human Little elongation complex subunit 2 protein (876-982AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Little elongation complex subunit 2 (Interactor of little elongator complex ELL subunit 2) (NMDA receptor-regulated protein 2), ICE2, BRCC1 NARG2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ICE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.