



PIK3R6 Antibody, Biotin conjugated

Product Code	CSB-PA719469LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5UE93
Immunogen	Recombinant Human Phosphoinositide 3-kinase regulatory subunit 6 protein (566-667AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoinositide 3-kinase regulatory subunit 6 (Phosphoinositide 3-kinase gamma adapter protein of 87 kDa) (p84 PI3K adapter protein) (p84 PIKAP) (p87 PI3K adapter protein) (p87PIKAP), PIK3R6, C17orf38
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIK3R6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.