



OBSCN Antibody, Biotin conjugated

Product Code	CSB-PA729422LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VST9
Immunogen	Recombinant Human Obscurin protein (1621-1712AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Obscurin (EC 2.7.11.1) (Obscurin-RhoGEF) (Obscurin-myosin light chain kinase) (Obscurin-MLCK), OBSCN, KIAA1556 KIAA1639
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OBSCN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.