



GPR107 Antibody, HRP conjugated

Product Code	CSB-PA736202LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VW38
Immunogen	Recombinant Human Protein GPR107 protein (151-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein GPR107 (Lung seven transmembrane receptor 1), GPR107, KIAA1624 LUSTR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR107
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.