



GPR107 Antibody, FITC conjugated

Product Code	CSB-PA736202LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VW38
Immunogen	Recombinant Human Protein GPR107 protein (151-262AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein GPR107 (Lung seven transmembrane receptor 1), GPR107, KIAA1624 LUSTR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR107