



UNC5D Antibody, HRP conjugated

Product Code	CSB-PA740925LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXZ4
Immunogen	Recombinant Human Netrin receptor UNC5D protein (90-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Netrin receptor UNC5D (Protein unc-5 homolog 4) (Protein unc-5 homolog D), UNC5D, KIAA1777 UNC5H4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UNC5D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.