




rbm24 Antibody

Product Code	CSB-PA745088LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76LC6
Immunogen	Recombinant Zebrafish RNA-binding protein 24 protein (2-78AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 24 (RNA-binding motif protein 24) (RNA-binding protein SEB-4), rbm24
Immunogen Species	Zebrafish
Research Area	Epigenetics and Nuclear Signaling
Target Names	rbm24

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: rbm24 antibody at 4.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41 kDa

Observed band size: 41 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.