



# rbm24 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA745088LC01DIL
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q76LC6
<b>Immunogen</b>	Recombinant Zebrafish RNA-binding protein 24 protein (2-78AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Zebrafish
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	RNA-binding protein 24 (RNA-binding motif protein 24) (RNA-binding protein SEB-4), rbm24
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	rbm24