



rbm24 Antibody, Biotin conjugated

Product Code	CSB-PA745088LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76LC6
Immunogen	Recombinant Zebrafish RNA-binding protein 24 protein (2-78AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 24 (RNA-binding motif protein 24) (RNA-binding protein SEB-4), rbm24
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Epigenetics and Nuclear Signaling
Target Names	rbm24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.