



CRTC3 Antibody, Biotin conjugated

Product Code	CSB-PA747623LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UUV7
Immunogen	Recombinant Human CREB-regulated transcription coactivator 3 protein (373-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CREB-regulated transcription coactivator 3 (Transducer of regulated cAMP response element-binding protein 3) (TORC-3) (Transducer of CREB protein 3), CRTC3, TORC3
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	CRTC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.