



# PNCK Antibody, Biotin conjugated

Product Code	CSB-PA754178LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P2M8
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase type 1B protein (1-115AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase type 1B (EC 2.7.11.17) (CaM kinase I beta) (CaM kinase IB) (CaM-KI beta) (CaMKI-beta) (Pregnancy up-regulated non-ubiquitously-expressed CaM kinase), PNCK
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PNCK