



ADGRG3 Antibody, HRP conjugated

Product Code	CSB-PA769807LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86Y34
Immunogen	Recombinant Human Adhesion G protein-coupled receptor G3 protein (98-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adhesion G protein-coupled receptor G3 (G-protein coupled receptor 97) (G-protein coupled receptor PGR26), ADGRG3, GPR97 PGR26
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADGRG3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.