



ADGRG3 Antibody, Biotin conjugated

Product Code	CSB-PA769807LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86Y34
Immunogen	Recombinant Human Adhesion G protein-coupled receptor G3 protein (98-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adhesion G protein-coupled receptor G3 (G-protein coupled receptor 97) (G-protein coupled receptor PGR26), ADGRG3, GPR97 PGR26
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADGRG3