



NTAN1 Antibody, Biotin conjugated

Product Code	CSB-PA822176LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AB6
Immunogen	Recombinant Human Protein N-terminal asparagine amidohydrolase protein (219-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein N-terminal asparagine amidohydrolase (EC 3.5.1.-) (Protein NH2-terminal asparagine amidohydrolase) (PNAA) (Protein NH2-terminal asparagine deamidase) (PNAD) (Protein N-terminal Asn amidase) (Protein N-terminal asparagine amidase) (Protein NTN-amidase), NTAN1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NTAN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.