



CELF6 Antibody, HRP conjugated

Product Code	CSB-PA822233LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96J87
Immunogen	Recombinant Human CUGBP Elav-like family member 6 protein (292-395AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CUGBP Elav-like family member 6 (CELF-6) (Bruno-like protein 6) (CUG-BP-and ETR-3-like factor 6) (RNA-binding protein BRUNOL-6), CELF6, BRUNOL6
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CELF6