



CASTOR3 Antibody, HRP conjugated

Product Code	CSB-PA822746LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NAP1
Immunogen	Recombinant Human Putative protein CASTOR 3 protein (1-163AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative protein GATS (STAG3 opposite strand transcript protein), GATS, STAG3OS
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GATS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.