



PIGT Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA836170LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q969N2 |
| Immunogen | Recombinant Human GPI transamidase component PIG-T protein (402-562AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | GPI transamidase component PIG-T (Phosphatidylinositol-glycan biosynthesis class T protein), PIGT |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | PIGT |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |