



PIGT Antibody, Biotin conjugated

Product Code	CSB-PA836170LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969N2
Immunogen	Recombinant Human GPI transamidase component PIG-T protein (402-562AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GPI transamidase component PIG-T (Phosphatidylinositol-glycan biosynthesis class T protein), PIGT
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PIGT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.