



FAR1 Antibody, Biotin conjugated

Product Code	CSB-PA837860LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVX9
Immunogen	Recombinant Human Fatty acyl-CoA reductase 1 protein (182-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty acyl-CoA reductase 1 (EC 1.2.1.84) (Male sterility domain-containing protein 2), FAR1, MLSTD2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	FAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.