



MIB2 Antibody, FITC conjugated

Product Code	CSB-PA839291LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AX9
Immunogen	Recombinant Human E3 ubiquitin-protein ligase MIB2 protein (1-128AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase MIB2 (EC 2.3.2.27) (Mind bomb homolog 2) (Novel zinc finger protein) (Novelzin) (Putative NF-kappa-B-activating protein 002N) (RING-type E3 ubiquitin transferase MIB2) (Skeletrophin) (Zinc finger ZZ type with ankyrin repeat domain protein 1), MIB2, SKD ZZANK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MIB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.