



MIB2 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA839291LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96AX9 |
| Immunogen | Recombinant Human E3 ubiquitin-protein ligase MIB2 protein (1-128AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | E3 ubiquitin-protein ligase MIB2 (EC 2.3.2.27) (Mind bomb homolog 2) (Novel zinc finger protein) (Novelzin) (Putative NF-kappa-B-activating protein 002N) (RING-type E3 ubiquitin transferase MIB2) (Skeletrophin) (Zinc finger ZZ type with ankyrin repeat domain protein 1), MIB2, SKD ZZANK1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | MIB2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |