



TMED6 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA840991LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8WW62 |
| Immunogen | Recombinant Human Transmembrane emp24 domain-containing protein 6 protein (22-170AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Transmembrane emp24 domain-containing protein 6 (p24 family protein gamma-5) (p24gamma5), TMED6 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | TMED6 |