



TMED6 Antibody, Biotin conjugated

Product Code	CSB-PA840991LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WW62
Immunogen	Recombinant Human Transmembrane emp24 domain-containing protein 6 protein (22-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane emp24 domain-containing protein 6 (p24 family protein gamma-5) (p24gamma5), TMED6
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TMED6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.