



SLC16A11 Antibody, Biotin conjugated

Product Code	CSB-PA843323LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCK7
Immunogen	Recombinant Human Monocarboxylate transporter 11 protein (428-471AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 11 (MCT 11) (Solute carrier family 16 member 11), SLC16A11, MCT11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC16A11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.