



RILP Antibody, Biotin conjugated

Product Code	CSB-PA850399LD01HU
Abbreviation	Rab-interacting lysosomal protein
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NA2
Immunogen	Recombinant Human Rab-interacting lysosomal protein (307-385AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Rab effector playing a role in late endocytic transport to degradative compartments. Involved in the regulation of lysosomal morphology and distribution. Induces recruitment of dynein-dynactin motor complexes to Rab7A-containing late endosome and lysosome compartments. Promotes centripetal migration of phagosomes and the fusion of phagosomes with the late endosomes and lysosomes.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab-interacting lysosomal protein, RILP
Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	RILP