



FCRL3 Antibody, FITC conjugated

Product Code	CSB-PA853475LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P31
Immunogen	Recombinant Human Fc receptor-like protein 3 protein (122-222AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fc receptor-like protein 3 (FcR-like protein 3) (FcRL3) (Fc receptor homolog 3) (FcRH3) (IFGP family protein 3) (hIFGP3) (Immune receptor translocation-associated protein 3) (SH2 domain-containing phosphatase anchor protein 2) (CD antigen CD307c), FCRL3, FCRH3 IFGP3 IRTA3 SPAP2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCRL3