



FCRL3 Antibody, Biotin conjugated

Product Code	CSB-PA853475LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P31
Immunogen	Recombinant Human Fc receptor-like protein 3 protein (122-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fc receptor-like protein 3 (FcR-like protein 3) (FcRL3) (Fc receptor homolog 3) (FcRH3) (IFGP family protein 3) (hIFGP3) (Immune receptor translocation-associated protein 3) (SH2 domain-containing phosphatase anchor protein 2) (CD antigen CD307c), FCRL3, FCRH3 IFGP3 IRTA3 SPAP2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCRL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.