



VN1R1 Antibody, HRP conjugated

Product Code	CSB-PA861973LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZP7
Immunogen	Recombinant Human Vomeronasal type-1 receptor 1 protein (1-56AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vomeronasal type-1 receptor 1 (G-protein coupled receptor GPCR24) (hGPCR24) (V1r-like receptor 1) (V3r-related gene protein) (Vomeronasal olfactory receptor chromosome 19 subtype I member 1), VN1R1, V1RL1 VNR1911
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VN1R1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.