



TRIM45 Antibody, Biotin conjugated

Product Code	CSB-PA862044LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8W5
Immunogen	Recombinant Human Tripartite motif-containing protein 45 protein (428-573AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 45 (RING finger protein 99), TRIM45, RNF99
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRIM45
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.