



MYB75 Antibody, HRP conjugated

Product Code	CSB-PA863725LB01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9FE25
Immunogen	Recombinant Arabidopsis thaliana Transcription factor MYB75 protein (91-167AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor MYB75 (Myb-related protein 75) (AtMYB75) (Production of anthocyanin pigment 1 protein) (Suc-induced anthocyanin accumulation locus 1), MYB75, PAP1 SIAA1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	MYB75
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.